

TECHNICAL SPECIFICATION
Cultivation Shelves

	CS 250/300_5_3	CS 250/200_7_4.2	CS 250/300_5_2.4	CS 250/200_7_3.4
LED light illumination	LED panel 100 x 60 cm		LED panel 80 x 60 cm	
External dimensions (W x D x H)	106 x 65 x 220 cm (height can be changed according to customer requirements)		86 x 65 x 220 cm (height can be changed according to customer requirements)	
Weight	75 kg	85 kg	75 kg	85 kg
Power	300 - 1200 W			
Warranty	1 year parts and labor			

Explanation to labeling: **VVV/XXX_Y_Z (CS 250/300_5_3)**

VVV - Maximum PPFD ($\mu\text{mol}\cdot\text{m}^{-2}\cdot\text{s}^{-1}$)

XXX - Distance of required PAR measurement [mm]

Y - Number of shelves

Z - Cultivation area (m^2)

Customized versions are also available.