

TECHNICAL SPECIFICATION

Walk-In Fytoscope FS-WI

LED light illumination	size on request
Maximum Illumination Intensity (Measured at 50 cm from the light panel)	Standard – 1,500 $\mu\text{mol}\cdot\text{m}^{-2}\cdot\text{s}^{-1}$ Optional – on request
Air flow	horizontal 0,25 $\text{m}\cdot\text{s}^{-1}$
Controlled Temperature Range	+10°C to +40°C (independent on light intensity; ambient room temperature up to +35 °C) optionally: - 5°C to 40°C
Humidification	ultrasonic humidifier (not hot steam)
Dehumidification	adsorption
Controlled Humidity Range	40 to 80 % (independent on light intensity)
External dimensions	depends on required space
Internal dimensions	depends on required space
Warranty	1 year parts and labor